

Quick Start Guide for PCR Tube Microcentrifuge

1. Supported Consumables

- **Primary Support:** This centrifuge is designed for **standard 1.5ml PCR centrifuge tubes**.
- **Compatibility:** Ensure your 1.5ml PCR tubes fit the rotor slots correctly. Do not use incompatible or deformed tubes.

1. Operating Steps

1. **Set Speed:** Turn the dial to set the desired speed (RPM).
2. **Set Duration:** Turn the dial to set the desired duration (seconds).
3. **Start:**
 - Press the "Start" button to begin centrifugation.
 - Wait for the run to finish.