

A Typical DNA Tailing Reaction

Overview

Protocol

1. Mix:

- 5.0 μ l (10X) TdT Buffer
- 5.0 μ l (2.5 mM) CoCl_2 solution provided
- 5.0 pmols DNA (330 ng for 100 bp, 1 μ g for 300 bp, 10 pmols DNA ends)*
- 0.5 μ l 10 mM dNTP (alpha- ^{32}P dATP may also be used)
- 0.5 μ l Terminal Transferase (20 units/ μ l)
- deionized water to a final volume of 50 μ l.

*To determine approximate amount of DNA (ng/pmol), multiply the number of base pairs by 0.66. Example: 300 bp x 0.66 = 198 ng/pmol. For 5.0 pmols multiply by 5, resulting in 990 ng/5 pmol.

The table below can be used as a guide (values are approximate and are given for a 30 minutes incubation at 37°C in the recommended buffer).

The rate of addition of dNTP's and thus the length of the tail is a function of the ratio of 3' DNA ends: dNTP concentration, and also which dNTP is used.

DNA Tailing Guide:

pmols 3' ends pmols dNTP	Tail Length			
	dA	dC	dG	dT
1:100	1-5	1-3	1-3	1-5
1:1,000	10-20	10-20	5-10	10-20
1:5,000	100-300	50-200	10-25	200-300

- Incubate at 37°C for 30 minutes.
- Stop the reaction by heating to 70°C for 10 minutes or by adding 10 μ l of 0.2 M EDTA (pH 8.0).